

CLEC 4E Polyclonal Antibody PE Conjugated

Catalog No: #C92697PE

Orders: order@signalwayantibody.comSupport: tech@signalwayantibody.com

Description

Product Name	CLEC 4E Polyclonal Antibody PE Conjugated
Host Species	Rabbit
Clonality	Polyclonal
Purification	affinity purified by Protein A
Applications	Flow-Cyt
Species Reactivity	Hu Ms Rat
Conjugates	PE
Target Name	CLEC 4E
Other Names	CLEC 4E
Accession No.	Uniprot:Q9ULY5Gene ID:26253
Uniprot	Q9ULY5
GeneID	26253
Calculated MW	25
Concentration	1 mg/ml
Formulation	0.01M TBS(pH7.4) with 1% BSA, 0.03% Proclin300 and 50% Glycerol.
Storage	Shipped at 4°C. Store at -20°C for one year. Avoid repeated freeze/thaw cycles.

Application Details

Flow-Cyt=2ug/Test

Note: This product is for in vitro research use only