

VGLUT3/SLC17A8 Polyclonal Antibody Biotin Conjugated

Catalog No: #C92717Bio

Orders: order@signalwayantibody.comSupport: tech@signalwayantibody.com

Description

Product Name	VGLUT3/SLC17A8 Polyclonal Antibody Biotin Conjugated
Host Species	Rabbit
Clonality	Polyclonal
Purification	affinity purified by Protein A
Applications	IHC ICC IF
Species Reactivity	Hu Ms Rat
Immunogen Type	1-100/589
Conjugates	Bio
Target Name	VGLUT3/SLC17A8
Other Names	VGLUT3/SLC17A8
Accession No.	Uniprot:Q8NDX2Gene ID:246213
Uniprot	Q8NDX2
GeneID	246213
Calculated MW	65
Formulation	Preservative: 15mM Sodium Azide, Constituents: 1% BSA, 0.01M PBS, pH 7.4
Storage	Shipped at 4°C. Store at -20°C for one year. Avoid repeated freeze/thaw cycles.

Application Details

IHC-P=1:50-200 IHC-F=1:50-200 ICC=1:50-200 IF=1:50-200

Note: This product is for in vitro research use only