

phospho-IRS1 (Tyr632) Polyclonal Antibody HRP Conjugated



Catalog No: #C92719HRP

Orders: order@signalwayantibody.comSupport: tech@signalwayantibody.com

Description

Product Name	phospho-IRS1 (Tyr632) Polyclonal Antibody HRP Conjugated
Host Species	Rabbit
Clonality	Polyclonal
Purification	affinity purified by Protein A
Applications	WB IHC-P IHC-F ICC
Species Reactivity	Hu Ms Rat
Immunogen Type	GD(p-Y)MP
Conjugates	HRP
Target Name	phospho-IRS1 (Tyr632)
Other Names	phospho-IRS1 (Tyr632)
Accession No.	Uniprot:P35568Gene ID:3667
Uniprot	P35568
GeneID	3667
Calculated MW	132
Formulation	Preservative: 15mM Sodium Azide, Constituents: 1% BSA, 0.01M PBS, pH 7.4
Storage	Shipped at 4°C. Store at -20°C for one year. Avoid repeated freeze/thaw cycles.

Application Details

WB=1:500-2000 IHC-P=1:50-200 IHC-F=1:50-200 ICC=1:50-200

Note: This product is for in vitro research use only