

# phospho-Smad2/Smad3 (Thr8) Polyclonal Antibody PE Conjugated



Catalog No: #C92723PE

Orders: [order@signalwayantibody.com](mailto:order@signalwayantibody.com)Support: [tech@signalwayantibody.com](mailto:tech@signalwayantibody.com)

## Description

Product Name	phospho-Smad2/Smad3 (Thr8) Polyclonal Antibody PE Conjugated
Host Species	Rabbit
Clonality	Polyclonal
Purification	affinity purified by Protein A
Applications	Flow-Cyt ICC IF
Species Reactivity	Hu Ms Rat
Immunogen Type	PF(p-T)PP
Conjugates	PE
Target Name	phospho-Smad2/Smad3 (Thr8)
Other Names	phospho-Smad2/Smad3 (Thr8)
Accession No.	Uniprot:P84022Gene ID:4087
Uniprot	P84022
GeneID	4087
Calculated MW	52
Formulation	0.01M TBS(pH7.4) with 1% BSA, 0.03% Proclin300 and 50% Glycerol.
Storage	Shipped at 4°C. Store at -20°C for one year. Avoid repeated freeze/thaw cycles.

## Application Details

Flow-Cyt=1:50-200 ICC=1:50-200 IF=1:50-200

Note: This product is for in vitro research use only