

Somatostatin 28 Polyclonal Antibody Biotin Conjugated

Catalog No: #C92731Bio

Orders: order@signalwayantibody.comSupport: tech@signalwayantibody.com

Description

Product Name	Somatostatin 28 Polyclonal Antibody Biotin Conjugated
Host Species	Rabbit
Clonality	Polyclonal
Purification	affinity purified by Protein A
Applications	IHC
Species Reactivity	Hu Ms Rat
Immunogen Type	51-116/116
Conjugates	Bio
Target Name	Somatostatin 28
Other Names	Somatostatin 28
Accession No.	Uniprot:P61278Gene ID:6750
Uniprot	P61278
GeneID	6750
Calculated MW	3
Formulation	0.01M TBS(pH7.4) with 1% BSA, 0.03% Proclin300 and 50% Glycerol.
Storage	Shipped at 4°C. Store at -20°C for one year. Avoid repeated freeze/thaw cycles.

Application Details

IHC=1:50-200

Note: This product is for in vitro research use only