

TRPM3 Polyclonal Antibody PE Conjugated

Catalog No: #C92738PE

Orders: order@signalwayantibody.comSupport: tech@signalwayantibody.com

Description

Product Name	TRPM3 Polyclonal Antibody PE Conjugated
Host Species	Rabbit
Clonality	Polyclonal
Purification	affinity purified by Protein A
Applications	IF
Species Reactivity	Hu Ms Rat
Immunogen Type	1121-1260/1732
Conjugates	PE
Target Name	TRPM3
Other Names	TRPM3
Accession No.	Uniprot:Q9HCF6Gene ID:80036
Uniprot	Q9HCF6
GeneID	80036
Calculated MW	197
Formulation	0.01M TBS(pH7.4) with 1% BSA, 0.03% Proclin300 and 50% Glycerol.
Storage	Shipped at 4°C. Store at -20°C for one year. Avoid repeated freeze/thaw cycles.

Application Details

IF=1:50-200

Note: This product is for in vitro research use only