

RNF8 Polyclonal Antibody APC Conjugated

Catalog No: #C92746APC

Orders: order@signalwayantibody.comSupport: tech@signalwayantibody.com

Description

Product Name	RNF8 Polyclonal Antibody APC Conjugated
Host Species	Rabbit
Clonality	Polyclonal
Purification	affinity purified by Protein A
Applications	IF
Species Reactivity	Hu Ms Rat
Immunogen Type	331-430/485
Conjugates	APC
Target Name	RNF8
Other Names	RNF8
Accession No.	Uniprot:O76064Gene ID:9025
Uniprot	O76064
GeneID	9025
Calculated MW	55
Formulation	0.01M TBS(pH7.4) with 1% BSA, 0.03% Proclin300 and 50% Glycerol.
Storage	Shipped at 4°C. Store at -20°C for one year. Avoid repeated freeze/thaw cycles.

Application Details

IF=1:50-200

Note: This product is for in vitro research use only