

## SCN5A/Nav1.5 Polyclonal Antibody PE Conjugated

Catalog No: #C92766PE

Orders: [order@signalwayantibody.com](mailto:order@signalwayantibody.com)Support: [tech@signalwayantibody.com](mailto:tech@signalwayantibody.com)

## Description

Product Name	SCN5A/Nav1.5 Polyclonal Antibody PE Conjugated
Host Species	Rabbit
Clonality	Polyclonal
Purification	affinity purified by Protein A
Applications	IF
Species Reactivity	Hu Ms Rat
Immunogen Type	51-150/2016
Conjugates	PE
Target Name	SCN5A/Nav1.5
Other Names	SCN5A/Nav1.5
Accession No.	Uniprot:Q14524Gene ID:6331
Uniprot	Q14524
GeneID	6331
Calculated MW	227
Formulation	0.01M TBS(pH7.4) with 1% BSA, 0.03% Proclin300 and 50% Glycerol.
Storage	Shipped at 4°C. Store at -20°C for one year. Avoid repeated freeze/thaw cycles.

## Application Details

IF=1:50-200

Note: This product is for in vitro research use only