

CTH Polyclonal Antibody Biotin Conjugated

Catalog No: #C92778Bio

Orders: order@signalwayantibody.comSupport: tech@signalwayantibody.com

Description

Product Name	CTH Polyclonal Antibody Biotin Conjugated
Host Species	Rabbit
Clonality	Polyclonal
Purification	affinity purified by Protein A
Applications	WB IHC IF
Species Reactivity	Hu Ms Rat
Immunogen Type	51-150/405
Conjugates	Bio
Target Name	CTH
Other Names	CTH
Accession No.	Uniprot:P32929Gene ID:1491
Uniprot	P32929
GeneID	1491
Calculated MW	45
Formulation	0.01M TBS(pH7.4) with 1% BSA, 0.03% Proclin300 and 50% Glycerol.
Storage	Shipped at 4°C. Store at -20°C for one year. Avoid repeated freeze/thaw cycles.

Application Details

WB=1:50-200 IHC-P=1:50-200 IHC-F=1:50-200 ICC=1:50-200 IF=1:50-200

Note: This product is for in vitro research use only