

STAT1 Polyclonal Antibody PE-Cy7 Conjugated

Catalog No: #C92784PE-Cy7

Orders: order@signalwayantibody.comSupport: tech@signalwayantibody.com

Description

Product Name	STAT1 Polyclonal Antibody PE-Cy7 Conjugated
Host Species	Rabbit
Clonality	Polyclonal
Purification	affinity purified by Protein A
Applications	ICC IF
Species Reactivity	Hu Ms Rat
Immunogen Type	251-350/750
Conjugates	PE-Cy7
Target Name	STAT1
Other Names	STAT1
Accession No.	Uniprot:P42224Gene ID:6772
Uniprot	P42224
GeneID	6772
Calculated MW	84
Formulation	0.01M TBS(pH7.4) with 1% BSA, 0.03% Proclin300 and 50% Glycerol.
Storage	Shipped at 4°C. Store at -20°C for one year. Avoid repeated freeze/thaw cycles.

Application Details

ICC=1:50-200 IF=1:50-200

Note: This product is for in vitro research use only