

TNF alpha(1F6) mouse Monoclonal Antibody Biotin Conjugated



Catalog No: #C92805Bio

Orders: order@signalwayantibody.comSupport: tech@signalwayantibody.com

Description

Product Name	TNF alpha(1F6) mouse Monoclonal Antibody Biotin Conjugated
Host Species	mouse
Clonality	Monoclonal
Purification	affinity purified by Protein G
Applications	WB IHC IF
Species Reactivity	Hu
Immunogen Type	77-233/233
Conjugates	Bio
Target Name	TNF alpha(1F6)
Other Names	TNF alpha(1F6)
Accession No.	Uniprot:P01375Gene ID:7124
Uniprot	P01375
GeneID	7124
Calculated MW	17/26
Formulation	0.01M TBS(pH7.4) with 1% BSA, 0.03% Proclin300 and 50% Glycerol.
Storage	Shipped at 4°C. Store at -20°C for one year. Avoid repeated freeze/thaw cycles.

Application Details

WB=1:50-200 IHC-P=1:50-200 IHC-F=1:50-200 ICC=1:50-200 IF=1:50-200

Note: This product is for in vitro research use only