

# IFN gamma(4E2) mouse Monoclonal Antibody Biotin Conjugated



Catalog No: #C92806Bio

Orders: [order@signalwayantibody.com](mailto:order@signalwayantibody.com)Support: [tech@signalwayantibody.com](mailto:tech@signalwayantibody.com)

## Description

Product Name	IFN gamma(4E2) mouse Monoclonal Antibody Biotin Conjugated
Host Species	mouse
Clonality	Monoclonal
Purification	affinity purified by Protein G
Applications	WB IHC IF
Species Reactivity	Hu
Immunogen Type	24-161/166 (C-6x His-Tag)
Conjugates	Bio
Target Name	IFN gamma(4E2)
Other Names	IFN gamma(4E2)
Accession No.	Uniprot:P01579Gene ID:3458
Uniprot	P01579
GeneID	3458
Calculated MW	17
Formulation	0.01M TBS(pH7.4) with 1% BSA, 0.03% Proclin300 and 50% Glycerol.
Storage	Shipped at 4°C. Store at -20°C for one year. Avoid repeated freeze/thaw cycles.

## Application Details

WB=1:50-200 IHC-P=1:50-200 IHC-F=1:50-200 ICC=1:50-200 IF=1:50-200

Note: This product is for in vitro research use only