

CRP(C11F2) mouse Monoclonal Antibody Gold Conjugated

Catalog No: #C92808Gold

Orders: order@signalwayantibody.comSupport: tech@signalwayantibody.com

Description

Product Name	CRP(C11F2) mouse Monoclonal Antibody Gold Conjugated
Host Species	mouse
Clonality	Monoclonal
Purification	affinity purified by Protein A
Applications	IEM IGS GICA
Species Reactivity	Hu
Conjugates	Gold
Target Name	CRP(C11F2)
Other Names	CRP(C11F2)
Accession No.	Uniprot:P02741Gene ID:1401
Uniprot	P02741
GeneID	1401
Calculated MW	25
Formulation	0.01M TBS(pH7.4) with 1% BSA, 0.03% Proclin300 and 50% Glycerol.
Storage	Shipped at 4°C. Store at -20°C for one year. Avoid repeated freeze/thaw cycles.

Application Details

IEM=1:20-200 IGS=1:20-200 GICA=1:20-200

Note: This product is for in vitro research use only