

PAG608 (3E4) mouse Monoclonal Antibody FITC Conjugated

Catalog No: #C92813FITC

Orders: order@signalwayantibody.comSupport: tech@signalwayantibody.com

Description

Product Name	PAG608 (3E4) mouse Monoclonal Antibody FITC Conjugated
Host Species	mouse
Clonality	Monoclonal
Purification	affinity purified by Protein A
Applications	IF
Species Reactivity	Hu Ms Rat
Immunogen Type	51-150/289
Conjugates	FITC
Target Name	PAG608 (3E4)
Other Names	PAG608 (3E4)
Accession No.	Uniprot:Q9HA38Gene ID:64393
Uniprot	Q9HA38
GeneID	64393
Calculated MW	32
Formulation	0.01M TBS(pH7.4) with 1% BSA, 0.03% Proclin300 and 50% Glycerol.
Storage	Shipped at 4°C. Store at -20°C for one year. Avoid repeated freeze/thaw cycles.

Application Details

IF=1:50-200

Note: This product is for in vitro research use only