

# HPV16 E6 protein mouse Monoclonal Antibody Biotin Conjugated

Catalog No: #C92815Bio

Orders: [order@signalwayantibody.com](mailto:order@signalwayantibody.com)Support: [tech@signalwayantibody.com](mailto:tech@signalwayantibody.com)

## Description

Product Name	HPV16 E6 protein mouse Monoclonal Antibody Biotin Conjugated
Host Species	mouse
Clonality	Monoclonal
Purification	affinity purified by Protein A
Applications	IHC IF
Species Reactivity	HPV16
Immunogen Type	85-158/158
Conjugates	Bio
Target Name	HPV16 E6 protein
Other Names	HPV16 E6 protein
Accession No.	Uniprot:P03126Gene ID:1489078
Uniprot	P03126
GeneID	1489078
Calculated MW	11
Formulation	0.01M TBS(pH7.4) with 1% BSA, 0.03% Proclin300 and 50% Glycerol.
Storage	Shipped at 4°C. Store at -20°C for one year. Avoid repeated freeze/thaw cycles.

## Application Details

IHC-P=1:50-200 IHC-F=1:50-200 IF=1:50-200

Note: This product is for in vitro research use only