

human IL-6 mouse Monoclonal Antibody HRP Conjugated

Catalog No: #C92819HRP

Orders: order@signalwayantibody.comSupport: tech@signalwayantibody.com

Description

Product Name	human IL-6 mouse Monoclonal Antibody HRP Conjugated
Host Species	mouse
Clonality	Monoclonal
Purification	affinity purified by Protein G
Applications	IHC ICC
Species Reactivity	Hu
Immunogen Type	30-212/212aa (C-6x His-Tag)
Conjugates	HRP
Target Name	human IL-6
Other Names	human IL-6
Accession No.	Uniprot:P05231Gene ID:3569
Uniprot	P05231
GeneID	3569
Calculated MW	23
Formulation	0.01M TBS(pH7.4) with 1% BSA, 0.03% Proclin300 and 50% Glycerol.
Storage	Shipped at 4°C. Store at -20°C for one year. Avoid repeated freeze/thaw cycles.

Application Details

IHC-P=1:50-200 IHC-F=1:50-200 ICC=1:50-200

Note: This product is for in vitro research use only