

G-CSF mouse Monoclonal Antibody HRP Conjugated

Catalog No: #C92820HRP

Orders: order@signalwayantibody.comSupport: tech@signalwayantibody.com

Description

Product Name	G-CSF mouse Monoclonal Antibody HRP Conjugated
Host Species	mouse
Clonality	Monoclonal
Purification	affinity purified by Protein A
Applications	WB IHC ICC
Species Reactivity	Hu Ms
Conjugates	HRP
Target Name	G-CSF
Other Names	G-CSF
Accession No.	Uniprot:P09919Gene ID:1440
Uniprot	P09919
GeneID	1440
Calculated MW	19
Formulation	0.01M TBS(pH7.4) with 1% BSA, 0.03% Proclin300 and 50% Glycerol.
Storage	Shipped at 4°C. Store at -20°C for one year. Avoid repeated freeze/thaw cycles.

Application Details

WB=1:500-2000 IHC-P=1:50-200 IHC-F=1:50-200 ICC=1:50-200

Note: This product is for in vitro research use only