

CD34 mouse Monoclonal Antibody PE Conjugated

Catalog No: #C92821PE

Orders: order@signalwayantibody.comSupport: tech@signalwayantibody.com

Description

Product Name	CD34 mouse Monoclonal Antibody PE Conjugated
Host Species	mouse
Clonality	Monoclonal
Purification	affinity purified by Protein G
Applications	ICC IF
Species Reactivity	Hu Rat
Immunogen Type	32-290/385
Conjugates	PE
Target Name	CD34
Other Names	CD34
Accession No.	Uniprot:P28906 Gene ID:947
Uniprot	P28906
GeneID	947
Calculated MW	39
Formulation	0.01M TBS(pH7.4) with 1% BSA, 0.03% Proclin300 and 50% Glycerol.
Storage	Shipped at 4°C. Store at -20°C for one year. Avoid repeated freeze/thaw cycles.

Application Details

ICC=1:50-200 IF=1:50-200

Note: This product is for in vitro research use only