

AFP(A2) mouse Monoclonal Antibody Biotin Conjugated

Catalog No: #C92823Bio

Orders: order@signalwayantibody.comSupport: tech@signalwayantibody.com

Description

Product Name	AFP(A2) mouse Monoclonal Antibody Biotin Conjugated
Host Species	mouse
Clonality	Monoclonal
Purification	affinity purified by Protein G
Applications	WB IHC IF
Species Reactivity	Hu Rat
Immunogen Type	full length
Conjugates	Bio
Target Name	AFP(A2)
Other Names	AFP(A2)
Accession No.	Uniprot:P02771Gene ID:174
Uniprot	P02771
GeneID	174
Calculated MW	65
Formulation	0.01M PBS(pH7.4) with 0.03% Proclin300.
Storage	Shipped at 4°C. Store at -20°C for one year. Avoid repeated freeze/thaw cycles.

Application Details

WB=1:50-200 IHC-P=1:50-200 IHC-F=1:50-200 ICC=1:50-200 IF=1:50-200

Note: This product is for in vitro research use only