

AFP(A4) mouse Monoclonal Antibody Biotin Conjugated

Catalog No: #C92824Bio

Orders: order@signalwayantibody.comSupport: tech@signalwayantibody.com

Description

| | |
|--------------------|---------------------------------------------------------------------------------|
| Product Name | AFP(A4) mouse Monoclonal Antibody Biotin Conjugated |
| Host Species | mouse |
| Clonality | Monoclonal |
| Purification | affinity purified by Protein A |
| Applications | WB IHC IF |
| Species Reactivity | Hu |
| Immunogen Type | full length |
| Conjugates | Bio |
| Target Name | AFP(A4) |
| Other Names | AFP(A4) |
| Accession No. | Uniprot:P02771Gene ID:174 |
| Uniprot | P02771 |
| GeneID | 174 |
| Calculated MW | 65 |
| Formulation | 0.01M PBS(pH7.4) with 0.03% Proclin300. |
| Storage | Shipped at 4°C. Store at -20°C for one year. Avoid repeated freeze/thaw cycles. |

Application Details

WB=1:50-200 IHC-P=1:50-200 IHC-F=1:50-200 ICC=1:50-200 IF=1:50-200

Note: This product is for in vitro research use only