

AFP(A4) mouse Monoclonal Antibody AP Conjugated

Catalog No: #C92824AP

Orders: order@signalwayantibody.comSupport: tech@signalwayantibody.com

Description

Product Name	AFP(A4) mouse Monoclonal Antibody AP Conjugated
Host Species	mouse
Clonality	Monoclonal
Purification	affinity purified by Protein A
Applications	WB IHC
Species Reactivity	Hu
Immunogen Type	full length
Conjugates	AP
Target Name	AFP(A4)
Other Names	AFP(A4)
Accession No.	Uniprot:P02771Gene ID:174
Uniprot	P02771
GeneID	174
Calculated MW	65
Formulation	0.01M PBS(pH7.4) with 0.03% Proclin300.
Storage	Shipped at 4°C. Store at -20°C for one year. Avoid repeated freeze/thaw cycles.

Application Details

WB=1:50-200 IHC-F=1:50-200

Note: This product is for in vitro research use only