

CEA(C3) mouse Monoclonal Antibody Biotin Conjugated

Catalog No: #C92825Bio

Orders: order@signalwayantibody.comSupport: tech@signalwayantibody.com

Description

Product Name	CEA(C3) mouse Monoclonal Antibody Biotin Conjugated
Host Species	mouse
Clonality	Monoclonal
Purification	affinity purified by Protein A
Applications	IHC IF
Species Reactivity	Hu
Immunogen Type	full length
Conjugates	Bio
Target Name	CEA(C3)
Other Names	CEA(C3)
Accession No.	Uniprot:P06731Gene ID:1048
Uniprot	P06731
GeneID	1048
Calculated MW	150-200
Formulation	0.01M PBS(pH7.4) with 0.03% Proclin300
Storage	Shipped at 4°C. Store at -20°C for one year. Avoid repeated freeze/thaw cycles.

Application Details

IHC-P=1:50-200 IHC-F=1:50-200 IF=1:50-200

Note: This product is for in vitro research use only