

# CEA(B5) mouse Monoclonal Antibody Gold Conjugated

Catalog No: #C92826Gold

Orders: [order@signalwayantibody.com](mailto:order@signalwayantibody.com)Support: [tech@signalwayantibody.com](mailto:tech@signalwayantibody.com)

## Description

Product Name	CEA(B5) mouse Monoclonal Antibody Gold Conjugated
Host Species	mouse
Clonality	Monoclonal
Purification	affinity purified by Protein A
Applications	IEM IGS GICA
Species Reactivity	Hu
Immunogen Type	full length
Conjugates	Gold
Target Name	CEA(B5)
Other Names	CEA(B5)
Accession No.	Uniprot:P06731Gene ID:1048
Uniprot	P06731
GeneID	1048
Calculated MW	150-200
Formulation	0.01M PBS(pH7.4) with 0.03% Proclin300
Storage	Shipped at 4°C. Store at -20°C for one year. Avoid repeated freeze/thaw cycles.

## Application Details

IEM=1:20-200 IGS=1:20-200 GICA=1:20-200

Note: This product is for in vitro research use only