

PSA (PA5) mouse Monoclonal Antibody Gold Conjugated

Catalog No: #C92828Gold

Orders: order@signalwayantibody.comSupport: tech@signalwayantibody.com

Description

Product Name	PSA (PA5) mouse Monoclonal Antibody Gold Conjugated
Host Species	mouse
Clonality	Monoclonal
Purification	affinity purified by Protein A
Applications	IEM IGS GICA
Species Reactivity	Hu
Immunogen Type	full length
Conjugates	Gold
Target Name	PSA (PA5)
Other Names	PSA (PA5)
Accession No.	Uniprot:P07288Gene ID:354
Uniprot	P07288
GeneID	354
Calculated MW	23
Formulation	0.01M PBS(pH7.4) with 0.03% Proclin300
Storage	Shipped at 4°C. Store at -20°C for one year. Avoid repeated freeze/thaw cycles.

Application Details

IEM=1:20-200 IGS=1:20-200 GICA=1:20-200

Note: This product is for in vitro research use only