

HBeAg (2E9) mouse Monoclonal Antibody Biotin Conjugated



Catalog No: #C92830Bio

Orders: order@signalwayantibody.comSupport: tech@signalwayantibody.com

Description

Product Name	HBeAg (2E9) mouse Monoclonal Antibody Biotin Conjugated
Host Species	mouse
Clonality	Monoclonal
Purification	affinity purified by Protein G
Applications	IHC
Species Reactivity	Hepatitis B virus
Conjugates	Bio
Target Name	HBeAg (2E9)
Other Names	HBeAg (2E9)
Accession No.	Uniprot:Q89714Gene ID:944568
Uniprot	Q89714
GeneID	944568
Calculated MW	18
Formulation	0.01M TBS(pH7.4) with 0.03% Proclin300
Storage	Shipped at 4°C. Store at -20°C for one year. Avoid repeated freeze/thaw cycles.

Application Details

IHC-F=1:50-200

Note: This product is for in vitro research use only