

HBeAg(2E6) mouse Monoclonal Antibody AP Conjugated

Catalog No: #C92831AP

Orders: order@signalwayantibody.comSupport: tech@signalwayantibody.com

Description

Product Name	HBeAg(2E6) mouse Monoclonal Antibody AP Conjugated
Host Species	mouse
Clonality	Monoclonal
Purification	affinity purified by Protein G
Applications	IHC
Species Reactivity	HBeAg
Conjugates	AP
Target Name	HBeAg(2E6)
Other Names	HBeAg(2E6)
Accession No.	Uniprot:Q89714Gene ID:944568
Uniprot	Q89714
GeneID	944568
Calculated MW	18
Formulation	0.01M TBS(pH7.4) with 0.03% Proclin300
Storage	Shipped at 4°C. Store at -20°C for one year. Avoid repeated freeze/thaw cycles.

Application Details

IHC-F=1:50-200

Note: This product is for in vitro research use only