

HCG Beta(4H7) mouse Monoclonal Antibody Biotin Conjugated

Catalog No: #C92833Bio

Orders: order@signalwayantibody.comSupport: tech@signalwayantibody.com

Description

Product Name	HCG Beta(4H7) mouse Monoclonal Antibody Biotin Conjugated
Host Species	mouse
Clonality	Monoclonal
Purification	affinity purified by Protein A
Applications	WB IHC IF
Species Reactivity	Hu
Conjugates	Bio
Target Name	HCG Beta(4H7)
Other Names	HCG Beta(4H7)
Accession No.	Uniprot:P01233Gene ID:1082
Uniprot	P01233
GeneID	1082
Calculated MW	16
Formulation	0.01M PBS(pH7.4) with 0.03% Proclin300
Storage	Shipped at 4°C. Store at -20°C for one year. Avoid repeated freeze/thaw cycles.

Application Details

WB=1:50-200 IHC-P=1:50-200 IHC-F=1:50-200 ICC=1:50-200 IF=1:50-200

Note: This product is for in vitro research use only