

HCG Beta(4H7) mouse Monoclonal Antibody Gold Conjugated



Catalog No: #C92833Gold

Orders: order@signalwayantibody.comSupport: tech@signalwayantibody.com

Description

Product Name	HCG Beta(4H7) mouse Monoclonal Antibody Gold Conjugated
Host Species	mouse
Clonality	Monoclonal
Purification	affinity purified by Protein A
Applications	IEM IGS GICA
Species Reactivity	Hu
Conjugates	Gold
Target Name	HCG Beta(4H7)
Other Names	HCG Beta(4H7)
Accession No.	Uniprot:P01233Gene ID:1082
Uniprot	P01233
GeneID	1082
Calculated MW	16
Formulation	0.01M PBS(pH7.4) with 0.03% Proclin300
Storage	Shipped at 4°C. Store at -20°C for one year. Avoid repeated freeze/thaw cycles.

Application Details

IEM=1:20-200 IGS=1:20-200 GICA=1:20-200

Note: This product is for in vitro research use only