

# MBP tag mouse Monoclonal Antibody HRP Conjugated

Catalog No: #C92842HRP

Orders: [order@signalwayantibody.com](mailto:order@signalwayantibody.com)Support: [tech@signalwayantibody.com](mailto:tech@signalwayantibody.com)

## Description

Product Name	MBP tag mouse Monoclonal Antibody HRP Conjugated
Host Species	mouse
Clonality	Monoclonal
Purification	affinity purified by Protein A
Applications	WB IHC ICC
Species Reactivity	MBP or MBP fusion proteins
Conjugates	HRP
Target Name	MBP tag
Other Names	MBP tag
Accession No.	Uniprot:P0AEX9Gene ID:1039564
Uniprot	P0AEX9
GeneID	1039564
Calculated MW	41
Formulation	0.01M TBS(pH7.4) with 1% BSA, 0.03% Proclin300 and 50% Glycerol.
Storage	Shipped at 4°C. Store at -20°C for one year. Avoid repeated freeze/thaw cycles.

## Application Details

WB=1:10000-20000 IHC-P=1:50-200 IHC-F=1:50-200 ICC=1:50-200

Note: This product is for in vitro research use only