

GFP-Tag mouse Monoclonal Antibody Biotin Conjugated

Catalog No: #C92846Bio

Orders: order@signalwayantibody.comSupport: tech@signalwayantibody.com

Description

Product Name	GFP-Tag mouse Monoclonal Antibody Biotin Conjugated
Host Species	mouse
Clonality	Monoclonal
Purification	affinity purified by Protein G
Applications	WB
Species Reactivity	Belt Jellyfish, GFP
Conjugates	Bio
Target Name	GFP-Tag
Other Names	GFP-Tag
Accession No.	Uniprot:P42212Gene ID:107331116
Uniprot	P42212
GeneID	107331116
Calculated MW	27
Formulation	0.01M TBS(pH7.4) with 1% BSA, 0.03% Proclin300 and 50% Glycerol.
Storage	Shipped at 4°C. Store at -20°C for one year. Avoid repeated freeze/thaw cycles.

Application Details

WB=1:50-200

Note: This product is for in vitro research use only