

Myc tag mouse Monoclonal Antibody HRP Conjugated

Catalog No: #C92847HRP

Orders: order@signalwayantibody.comSupport: tech@signalwayantibody.com

Description

| | |
|--------------------|---|
| Product Name | Myc tag mouse Monoclonal Antibody HRP Conjugated |
| Host Species | mouse |
| Clonality | Monoclonal |
| Purification | affinity purified by Protein G |
| Applications | WB ICC |
| Species Reactivity | Hu Myc tag |
| Conjugates | HRP |
| Target Name | Myc tag |
| Other Names | Myc tag |
| Accession No. | Uniprot:P01106Gene ID:4609 |
| Uniprot | P01106 |
| GeneID | 4609 |
| Formulation | 0.01M TBS(pH7.4) with 1% BSA, 0.03% Proclin300 and 50% Glycerol. |
| Storage | Shipped at 4°C. Store at -20°C for one year. Avoid repeated freeze/thaw cycles. |

Application Details

WB=1:10000-20000 ICC=1:50-200

Note: This product is for in vitro research use only