

Myc tag mouse Monoclonal Antibody AP Conjugated

Catalog No: #C92847AP

Orders: order@signalwayantibody.comSupport: tech@signalwayantibody.com

Description

Product Name	Myc tag mouse Monoclonal Antibody AP Conjugated
Host Species	mouse
Clonality	Monoclonal
Purification	affinity purified by Protein G
Applications	WB ICC
Species Reactivity	Hu Myc tag
Conjugates	AP
Target Name	Myc tag
Other Names	Myc tag
Accession No.	Uniprot:P01106Gene ID:4609
Uniprot	P01106
GeneID	4609
Formulation	0.01M TBS(pH7.4) with 1% BSA, 0.03% Proclin300 and 50% Glycerol.
Storage	Shipped at 4°C. Store at -20°C for one year. Avoid repeated freeze/thaw cycles.

Application Details

WB=1:10000-20000 ICC=1:50-200

Note: This product is for in vitro research use only