

## Myc tag mouse Monoclonal Antibody PE Conjugated

Catalog No: #C92847PE

Orders: [order@signalwayantibody.com](mailto:order@signalwayantibody.com)Support: [tech@signalwayantibody.com](mailto:tech@signalwayantibody.com)

## Description

Product Name	Myc tag mouse Monoclonal Antibody PE Conjugated
Host Species	mouse
Clonality	Monoclonal
Purification	affinity purified by Protein G
Applications	Flow-Cyt ICC IF
Species Reactivity	Hu Myc tag
Conjugates	PE
Target Name	Myc tag
Other Names	Myc tag
Accession No.	Uniprot:P01106Gene ID:4609
Uniprot	P01106
GeneID	4609
Formulation	0.01M TBS(pH7.4) with 1% BSA, 0.03% Proclin300 and 50% Glycerol.
Storage	Shipped at 4°C. Store at -20°C for one year. Avoid repeated freeze/thaw cycles.

## Application Details

Flow-Cyt=2ug/Test ICC=1:50-200 IF=1:50-200

Note: This product is for in vitro research use only