

alpha Tubulin mouse Monoclonal Antibody FITC Conjugated



Catalog No: #C92851FITC

Orders: order@signalwayantibody.comSupport: tech@signalwayantibody.com

Description

Product Name	alpha Tubulin mouse Monoclonal Antibody FITC Conjugated
Host Species	mouse
Clonality	Monoclonal
Purification	affinity purified by Protein G
Applications	ICC IF
Species Reactivity	Hu Ms Rat
Conjugates	FITC
Target Name	alpha Tubulin
Other Names	alpha Tubulin
Accession No.	Uniprot:P68366Gene ID:7277
Uniprot	P68366
GeneID	7277
Calculated MW	50
Formulation	0.01M TBS(pH7.4) with 1% BSA, 0.03% Proclin300 and 50% Glycerol.
Storage	Shipped at 4°C. Store at -20°C for one year. Avoid repeated freeze/thaw cycles.

Application Details

ICC=1:50-200 IF=1:50-200

Note: This product is for in vitro research use only