

beta I Tubulin(Loading Control) mouse Monoclonal Antibody PE-Cy5 Conjugated



Catalog No: #C92852PE-Cy5

Orders: order@signalwayantibody.comSupport: tech@signalwayantibody.com

Description

Product Name	beta I Tubulin(Loading Control) mouse Monoclonal Antibody PE-Cy5 Conjugated
Host Species	mouse
Clonality	Monoclonal
Purification	affinity purified by Protein G
Applications	Flow-Cyt ICC IF
Species Reactivity	Hu Ms Rat
Conjugates	PE-Cy5
Target Name	beta I Tubulin(Loading Control)
Other Names	beta I Tubulin(Loading Control)
Accession No.	Uniprot:Q9H4B7Gene ID:81027
Uniprot	Q9H4B7
GeneID	81027
Calculated MW	50
Formulation	0.01M TBS(pH7.4) with 1% BSA, 0.03% Proclin300 and 50% Glycerol.
Storage	Shipped at 4°C. Store at -20°C for one year. Avoid repeated freeze/thaw cycles.

Application Details

Flow-Cyt=1:50-200 ICC=1:50-200 IF=1:50-200

Note: This product is for in vitro research use only