

## EGFR mouse Monoclonal Antibody Gold Conjugated

Catalog No: #C92855Gold

Orders: [order@signalwayantibody.com](mailto:order@signalwayantibody.com)Support: [tech@signalwayantibody.com](mailto:tech@signalwayantibody.com)

## Description

Product Name	EGFR mouse Monoclonal Antibody Gold Conjugated
Host Species	mouse
Clonality	Monoclonal
Purification	affinity purified by Protein G
Applications	WB ICC IF
Species Reactivity	Hu Ms Rat
Conjugates	Gold
Target Name	EGFR
Other Names	EGFR
Accession No.	Uniprot:P00533Gene ID:1956
Uniprot	P00533
GeneID	1956
Calculated MW	130
Formulation	0.01M TBS(pH7.4) with 1% BSA, 0.03% Proclin300 and 50% Glycerol.
Storage	Shipped at 4°C. Store at -20°C for one year. Avoid repeated freeze/thaw cycles.

## Application Details

WB=1:50-200 ICC=1:50-200 IF=1:50-200

Note: This product is for in vitro research use only