

EGFR mouse Monoclonal Antibody FITC Conjugated

Catalog No: #C92855FITC

Orders: order@signalwayantibody.comSupport: tech@signalwayantibody.com

Description

| | |
|--------------------|---|
| Product Name | EGFR mouse Monoclonal Antibody FITC Conjugated |
| Host Species | mouse |
| Clonality | Monoclonal |
| Purification | affinity purified by Protein G |
| Applications | Flow-Cyt IF |
| Species Reactivity | Hu Ms Rat |
| Conjugates | FITC |
| Target Name | EGFR |
| Other Names | EGFR |
| Accession No. | Uniprot:P00533Gene ID:1956 |
| Uniprot | P00533 |
| GeneID | 1956 |
| Calculated MW | 130 |
| Formulation | 0.01M TBS(pH7.4) with 1% BSA, 0.03% Proclin300 and 50% Glycerol. |
| Storage | Shipped at 4°C. Store at -20°C for one year. Avoid repeated freeze/thaw cycles. |

Application Details

Flow-Cyt=2ug/Test IF=1:100-500

Note: This product is for in vitro research use only