

C-erbB-2/HER2 mouse Monoclonal Antibody Cy3 Conjugated

Catalog No: #C92856Cy3

Orders: order@signalwayantibody.comSupport: tech@signalwayantibody.com

Description

Product Name	C-erbB-2/HER2 mouse Monoclonal Antibody Cy3 Conjugated
Host Species	mouse
Clonality	Monoclonal
Purification	affinity purified by Protein G
Applications	IF
Species Reactivity	Hu Ms Rat
Immunogen Type	full length
Conjugates	Cy3
Target Name	C-erbB-2/HER2
Other Names	C-erbB-2/HER2
Accession No.	Uniprot:P04626Gene ID:2064
Uniprot	P04626
GeneID	2064
Calculated MW	138
Formulation	0.01M TBS(pH7.4) with 1% BSA, 0.03% Proclin300 and 50% Glycerol.
Storage	Shipped at 4°C. Store at -20°C for one year. Avoid repeated freeze/thaw cycles.

Application Details

IF=1:50-200

Note: This product is for in vitro research use only