

## Ki-67 mouse Monoclonal Antibody FITC Conjugated

Catalog No: #C92863FITC

Orders: [order@signalwayantibody.com](mailto:order@signalwayantibody.com)Support: [tech@signalwayantibody.com](mailto:tech@signalwayantibody.com)

## Description

Product Name	Ki-67 mouse Monoclonal Antibody FITC Conjugated
Host Species	mouse
Clonality	Monoclonal
Purification	affinity purified by Protein G
Applications	Flow-Cyt ICC IF
Species Reactivity	Hu Ms Rat
Conjugates	FITC
Target Name	Ki-67
Other Names	Ki-67
Accession No.	Uniprot:P46013Gene ID:4288
Uniprot	P46013
GeneID	4288
Calculated MW	358
Formulation	0.01M TBS(pH7.4) with 1% BSA, 0.03% Proclin300 and 50% Glycerol.
Storage	Shipped at 4°C. Store at -20°C for one year. Avoid repeated freeze/thaw cycles.

## Application Details

Flow-Cyt=2ug/Test ICC=1:50-200 IF=1:50-200

Note: This product is for in vitro research use only