

## NSE mouse Monoclonal Antibody HRP Conjugated

Catalog No: #C92864HRP

Orders: [order@signalwayantibody.com](mailto:order@signalwayantibody.com)Support: [tech@signalwayantibody.com](mailto:tech@signalwayantibody.com)

## Description

Product Name	NSE mouse Monoclonal Antibody HRP Conjugated
Host Species	mouse
Clonality	Monoclonal
Purification	affinity purified by Protein G
Applications	WB
Species Reactivity	Hu Ms Rat
Conjugates	HRP
Target Name	NSE
Other Names	NSE
Accession No.	Uniprot:P09104Gene ID:2026
Uniprot	P09104
GeneID	2026
Calculated MW	48
Formulation	0.01M TBS(pH7.4) with 1% BSA, 0.03% Proclin300 and 50% Glycerol.
Storage	Shipped at 4°C. Store at -20°C for one year. Avoid repeated freeze/thaw cycles.

## Application Details

WB=1:500-2000

Note: This product is for in vitro research use only