

AFP mouse Monoclonal Antibody Biotin Conjugated

Catalog No: #C92868Bio

Orders: order@signalwayantibody.comSupport: tech@signalwayantibody.com

Description

Product Name	AFP mouse Monoclonal Antibody Biotin Conjugated
Host Species	mouse
Clonality	Monoclonal
Purification	affinity purified by Protein G
Applications	WB IHC ICC IF
Species Reactivity	Hu
Conjugates	Bio
Target Name	AFP
Other Names	AFP
Accession No.	Uniprot:P02771Gene ID:174
Uniprot	P02771
GeneID	174
Calculated MW	70
Formulation	0.01M TBS(pH7.4) with 1% BSA, 0.03% Proclin300 and 50% Glycerol.
Storage	Shipped at 4°C. Store at -20°C for one year. Avoid repeated freeze/thaw cycles.

Application Details

WB=1:50-200 IHC-F=1:50-200 ICC=1:50-200 IF=1:50-200

Note: This product is for in vitro research use only