

Progesterone Receptor mouse Monoclonal Antibody PE Conjugated



Catalog No: #C92879PE

Orders: order@signalwayantibody.comSupport: tech@signalwayantibody.com

Description

Product Name	Progesterone Receptor mouse Monoclonal Antibody PE Conjugated
Host Species	mouse
Clonality	Monoclonal
Purification	affinity purified by Protein G
Applications	ICC IF
Species Reactivity	Hu Ms Rat
Conjugates	PE
Target Name	Progesterone Receptor
Other Names	Progesterone Receptor
Accession No.	Uniprot:P06401Gene ID:5241
Uniprot	P06401
GeneID	5241
Calculated MW	103
Formulation	0.01M TBS(pH7.4) with 1% BSA, 0.03% Proclin300 and 50% Glycerol.
Storage	Shipped at 4°C. Store at -20°C for one year. Avoid repeated freeze/thaw cycles.

Application Details

ICC=1:50-200 IF=1:50-200

Note: This product is for in vitro research use only