

TNF alpha mouse Monoclonal Antibody PE-Cy3 Conjugated

Catalog No: #C92887PE-Cy3

Orders: order@signalwayantibody.comSupport: tech@signalwayantibody.com

Description

Product Name	TNF alpha mouse Monoclonal Antibody PE-Cy3 Conjugated
Host Species	mouse
Clonality	Monoclonal
Purification	affinity purified by Protein G
Applications	ICC IF
Species Reactivity	Hu Ms Rat
Conjugates	PE-Cy3
Target Name	TNF alpha
Other Names	TNF alpha
Accession No.	Uniprot:P01375Gene ID:7124
Uniprot	P01375
GeneID	7124
Calculated MW	17/26
Formulation	0.01M TBS(pH7.4) with 1% BSA, 0.03% Proclin300 and 50% Glycerol.
Storage	Shipped at 4°C. Store at -20°C for one year. Avoid repeated freeze/thaw cycles.

Application Details

ICC=1:50-200 IF=1:50-200

Note: This product is for in vitro research use only