

# TNF alpha mouse Monoclonal Antibody HRP Conjugated

Catalog No: #C92888HRP

Orders: [order@signalwayantibody.com](mailto:order@signalwayantibody.com)Support: [tech@signalwayantibody.com](mailto:tech@signalwayantibody.com)

## Description

Product Name	TNF alpha mouse Monoclonal Antibody HRP Conjugated
Host Species	mouse
Clonality	Monoclonal
Purification	affinity purified by Protein G
Applications	WB IHC ICC
Species Reactivity	Hu Ms Rat
Conjugates	HRP
Target Name	TNF alpha
Other Names	TNF alpha
Accession No.	Uniprot:P01375Gene ID:7124
Uniprot	P01375
GeneID	7124
Calculated MW	17/26
Formulation	0.01M TBS(pH7.4) with 1% BSA, 0.03% Proclin300 and 50% Glycerol.
Storage	Shipped at 4°C. Store at -20°C for one year. Avoid repeated freeze/thaw cycles.

## Application Details

WB=1:10000-20000 IHC-P=1:50-200 IHC-F=1:50-200 ICC=1:50-200

Note: This product is for in vitro research use only