

Transferrin(Serum Loading Control) mouse Monoclonal Antibody HRP Conjugated

Catalog No: #C92889HRP

Orders: order@signalwayantibody.comSupport: tech@signalwayantibody.com

Description

Product Name	Transferrin(Serum Loading Control) mouse Monoclonal Antibody HRP Conjugated
Host Species	mouse
Clonality	Monoclonal
Purification	affinity purified by Protein G
Applications	WB IHC ICC
Species Reactivity	Hu
Conjugates	HRP
Target Name	Transferrin(Serum Loading Control)
Other Names	Transferrin(Serum Loading Control)
Accession No.	Uniprot:P02787Gene ID:7018
Uniprot	P02787
GeneID	7018
Calculated MW	77
Formulation	0.01M TBS(pH7.4) with 1% BSA, 0.03% Proclin300 and 50% Glycerol.
Storage	Shipped at 4°C. Store at -20°C for one year. Avoid repeated freeze/thaw cycles.

Application Details

WB=1:10000-20000 IHC-P=1:50-200 IHC-F=1:50-200 ICC=1:50-200

Note: This product is for in vitro research use only