

## SGK1 mouse Monoclonal Antibody HRP Conjugated

Catalog No: #C92893HRP

Orders: [order@signalwayantibody.com](mailto:order@signalwayantibody.com)Support: [tech@signalwayantibody.com](mailto:tech@signalwayantibody.com)

## Description

Product Name	SGK1 mouse Monoclonal Antibody HRP Conjugated
Host Species	mouse
Clonality	Monoclonal
Purification	affinity purified by Protein G
Applications	IHC ICC
Species Reactivity	Hu Ms Rat
Conjugates	HRP
Target Name	SGK1
Other Names	SGK1
Accession No.	Uniprot:O00141Gene ID:6446
Uniprot	O00141
GeneID	6446
Calculated MW	49
Formulation	0.01M TBS(pH7.4) with 1% BSA, 0.03% Proclin300 and 50% Glycerol.
Storage	Shipped at 4°C. Store at -20°C for one year. Avoid repeated freeze/thaw cycles.

## Application Details

IHC-P=1:50-200 IHC-F=1:50-200 ICC=1:50-200

Note: This product is for in vitro research use only