

SGK1 mouse Monoclonal Antibody PE Conjugated

Catalog No: #C92893PE

Orders: order@signalwayantibody.comSupport: tech@signalwayantibody.com

Description

Product Name	SGK1 mouse Monoclonal Antibody PE Conjugated
Host Species	mouse
Clonality	Monoclonal
Purification	affinity purified by Protein G
Applications	ICC IF
Species Reactivity	Hu Ms Rat
Conjugates	PE
Target Name	SGK1
Other Names	SGK1
Accession No.	Uniprot:O00141Gene ID:6446
Uniprot	O00141
GeneID	6446
Calculated MW	49
Formulation	0.01M TBS(pH7.4) with 1% BSA, 0.03% Proclin300 and 50% Glycerol.
Storage	Shipped at 4°C. Store at -20°C for one year. Avoid repeated freeze/thaw cycles.

Application Details

ICC=1:50-200 IF=1:50-200

Note: This product is for in vitro research use only