

## Bax mouse Monoclonal Antibody PE-Cy7 Conjugated

Catalog No: #C92895PE-Cy7

Orders: [order@signalwayantibody.com](mailto:order@signalwayantibody.com)Support: [tech@signalwayantibody.com](mailto:tech@signalwayantibody.com)

## Description

Product Name	Bax mouse Monoclonal Antibody PE-Cy7 Conjugated
Host Species	mouse
Clonality	Monoclonal
Purification	affinity purified by Protein G
Applications	ICC IF
Species Reactivity	Hu Ms Rat
Conjugates	PE-Cy7
Target Name	Bax
Other Names	Bax
Accession No.	Uniprot:Q07812Gene ID:581
Uniprot	Q07812
GeneID	581
Calculated MW	21
Formulation	0.01M TBS(pH7.4) with 1% BSA, 0.03% Proclin300 and 50% Glycerol.
Storage	Shipped at 4°C. Store at -20°C for one year. Avoid repeated freeze/thaw cycles.

## Application Details

ICC=1:50-200 IF=1:50-200

Note: This product is for in vitro research use only