

Caspase-3 mouse Monoclonal Antibody Cy3 Conjugated

Catalog No: #C92897Cy3

Orders: order@signalwayantibody.comSupport: tech@signalwayantibody.com

Description

Product Name	Caspase-3 mouse Monoclonal Antibody Cy3 Conjugated
Host Species	mouse
Clonality	Monoclonal
Purification	affinity purified by Protein G
Applications	ICC IF
Species Reactivity	Hu Ms Rat
Immunogen Type	1-100/277
Conjugates	Cy3
Target Name	Caspase-3
Other Names	Caspase-3
Accession No.	Uniprot:P42574Gene ID:836
Uniprot	P42574
GeneID	836
Calculated MW	17/28
Formulation	0.01M TBS(pH7.4) with 1% BSA, 0.03% Proclin300 and 50% Glycerol.
Storage	Shipped at 4°C. Store at -20°C for one year. Avoid repeated freeze/thaw cycles.

Application Details

ICC=1:50-200 IF=1:50-200

Note: This product is for in vitro research use only