

GSK-3 Beta mouse Monoclonal Antibody HRP Conjugated



Catalog No: #C92898HRP

Orders: order@signalwayantibody.comSupport: tech@signalwayantibody.com

Description

| | |
|--------------------|---|
| Product Name | GSK-3 Beta mouse Monoclonal Antibody HRP Conjugated |
| Host Species | mouse |
| Clonality | Monoclonal |
| Purification | affinity purified by Protein G |
| Applications | WB IHC ICC |
| Species Reactivity | Hu Ms Rat |
| Immunogen Type | 1-100/420 |
| Conjugates | HRP |
| Target Name | GSK-3 Beta |
| Other Names | GSK-3 Beta |
| Accession No. | Uniprot:P49841Gene ID:2932 |
| Uniprot | P49841 |
| GeneID | 2932 |
| Calculated MW | 47 |
| Formulation | 0.01M TBS(pH7.4) with 1% BSA, 0.03% Proclin300 and 50% Glycerol. |
| Storage | Shipped at 4°C. Store at -20°C for one year. Avoid repeated freeze/thaw cycles. |

Application Details

WB=1:10000-20000 IHC-P=1:50-200 IHC-F=1:50-200 ICC=1:50-200

Note: This product is for in vitro research use only